

Diffusion Weighted Imaging (DWI)

– ADC

BRAIN

[Biomarker](#) / Brain / Diffusion Weighted Imaging (DWI) – ADC

DESCRIPTION

Quibim’s app enables to automatically calculate the apparent diffusion coefficient (ADC) to quantify the diffusion of water molecules in tissues.

The app provides information on the cellular and microstructural organization of tissues and highlights the relationship between pathological changes and the diffusion of water molecules.

This app is useful for the characterization of tumoral processes to improve early detection, diagnosis, treatment response evaluation and follow-up of patients with gliomas, meningiomas and lymphomas as well as ischemic injuries.

DISEASES

Brain cancer including:

- Gliomas
- Meningiomas
- Lymphomas
- Ischemic injuries

REFERENCES



[Martí-Bonmatí L, Ramírez-Fuentes C, Alberich-Bayarri Á, Ruiz-Llorca C. State-of-the-art of bone marrow imaging in multiple myeloma. Curr Opin Oncol. 2015 Nov; 27\(6\): 540-50.](#)

SHARE:   





Imaging Center	PAT5	Patient Name	Glioblastoma
Modality	MR	Patient ID	

Study Description	RM111 - RM CEREBRAL	Patient Sex	F
Study Date	11/04/2016	Birthdate	



Min-Max b-values



	ADC [mm ² /s]
Mean ± Std	1.22e-03 ± 1.75e-03
Median	9.95e-04
P25	8.19e-04
P75	1.42e-03

QUIBIM S.L. - Quantitative Imaging Biomarkers in Medicine
EDIFICIO EUROPA - Avenida Aragon 30, 13th Floor, Office J. Valencia (SPAIN)
contact@quibim.com | www.quibim.com | Phone: +34 961 243 225

